Stat Profile®/Nova Hematocrit Controls

CONTROL 1 2

LOT

| | | |
|-----------|--------|---------|
| CONTROL 1 | 902400 | 2012-02 |
| CONTROL 2 | 902401 | 2012-02 |

| | | Expected Ranges, Zu erwartende Bereiche, Αναμενόμενη περιοχή τιμών, Gamas esperadas, Plages prévues, Intervalli attesi, Intervalo Esperado, Förväntade områden | |
|--------------------------|---|--|-----------|
| | | CONTROL 1 | CONTROL 2 |
| Stat Profile® Hematocrit | % | 31-35 | 50-54 |
| Nova Hematocrit | % | 31-35 | 50-54 |

Product Description: An aqueous quality control material for monitoring the performance of hematocrit. For use on Nova Biomedical analyzers only. The controls are formulated at two levels:

CONTROL 1 Control 1: Normal, Low Values
CONTROL 2 Control 2: Normal, High Values

Intended Use: For in vitro diagnostic use for monitoring the performance of hematocrit on Stat Profile and Nova Analyzers.

Methodology: Refer to Stat Profile or Nova Instructions for Use Manual for Methodology and Principals.

Composition: An electrolyte solution. The concentrations of the active ingredients have been adjusted to give a conductivity signal equivalent to a known hematocrit value in whole blood and to take into account the effect of sodium concentration on the conductivity on the samples. Each ampoule contains 1.7 mL. Contains no constituents of human origin, however, good laboratory practices should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2)

Warnings and Cautions: Intended for In Vitro Diagnostic use. Refer to Stat Profile or Nova Instructions for Use Manual for complete directions. DO NOT FREEZE. Follow standard practices for handling laboratory reagents.

Storage: Store at 15°C - 30°C. DO NOT FREEZE. An expiration date is printed on each ampoule.

Directions for Use: Verify that the lot number appearing on the expected ranges table corresponds to the lot on the ampoule. Snap open the ampoule (protecting fingers with gauze or gloves.) Aspirate liquid from the ampoule to the analyzer. Refer to Stat Profile or Nova Instructions for Use Manual for complete directions.

Limitations: Hematocrit controls are intended for monitoring the performance of Nova Biomedical analyzers ONLY. Not for use on any other manufacturers' analyzers.

Traceability of Standards: Hematocrit values are traceable to NIST Standards.

Expected Ranges: The expected range for each analyte was determined at Nova Biomedical by using multiple runs of each level of control on multiple analyzers. The expected range indicates the maximum deviations from the mean value which may be expected under differing laboratory conditions for analyzers operating within specifications. Refer to expected ranges table.

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Waltham, MA 02454 U.S.A.